



# **Small Headless Screw**

Surgical Technique | TriMed Cannulated Screw System








Screws and their respective instrumentation are color coded by screw diameter. See page 4 for size and color reference chart.

# Wire/Drill Guide Assembly


See technique on page 3 for 1.7mm screws

- Snap Modular Handle into Drill Cannula.
- Slide Wire Guide into Drill Cannula until fully seated.




# K-wire Insertion

- Drive the appropriate size K-wire through the guide to desired depth.
- Measure K-wire depth through the guide window. (See technique on page 3 for 3.5mm screws)
- If desired, advance K-wire further to help prevent disengagement when drilling over K-wire.



# **Wire Guide Removal**


- Withdraw the Wire Guide from the Drill Cannula.
- Select the corresponding drill size for the intended screw diameter.



# **Site Preparation**

- Drill to the desired depth over the K-wire.
- The depth of the hole can be checked through the guide window.
- · Remove the drill bit and Drill Cannula.
- Countersink hole as needed to recess the screw head within the cortical bone.





# **Screw Insertion**

- Select the appropriate screw length.
- Drive screw to desired position and remove K-wire.



Indications, contraindications, warnings and precautions related to TriMed Compression Screws reference IFU on trimedortho.com/ifu

#### **TECHNIQUE**


# 1.7mm Screws:



Note: This size screw does not require a drill and does not use a drill cannula.

- Snap Modular Handle Into the Wire Guide.
- Drive K-wire through the guide to desired depth.
- Measure K-wire depth through the guide window.
- Remove Wire Guide from K-wire.
- Insert screw (as illustrated in step 5).

# 3.5mm Screws:




# **Upper Mark:**

For lengths 35, 40 and 45mm. The screw size is determined by adding **20mm** to upper mark measurement.

### **Lower Mark:**

For lengths 34mm or less.



# All implants made from surgical grade titanium













Screw	Length	Thread	Head	Wire Guide	Drill Cannula	K-wire	Drill Bit	Countersink
<b>1.7</b> L17xx	08–14mm <sup>1</sup>	1.7mm	2.4mm	WGUIDE-1.7	n/a	WIRE-0.7/080	[ self-drilling ]	HSINK-1.7
<b>2.3</b> L23xx	10–20mm² 20–26mm¹ 26-28mm²	2.3mm	3.0mm	WGUIDE-2.3	CANNULA-2.3	WIRE-0.8/120	DRILL-1.6/095C	HSINK-2.3
<b>3.0</b> L30xx	10–20mm² 20–26mm¹ 26-36mm²	3.0mm	4.0mm	WGUIDE-3.0	CANNULA-3.0	WIRE-1.1/120	DRILL-2.1/110C	HSINK-3.0
<b>3.5</b> L35xx	20–32mm² 35–45mm⁵	3.5mm	4.5mm	WGUIDE-3.5	CANNULA-3.5	WIRE-1.1/120	DRILL-2.4/120C	HSINK-3.5

 $mm^1 = 1mm$  increments

mm<sup>2</sup> = 2mm increments

 $mm^5 = 5mm$  increments



TriMed, Inc. / 27533 Avenue Hopkins / Valencia, CA 91355 USA / 800-633-7221 / www.trimedortho.com

The technique presented is one suggested surgical technique. The decision to use a specific implant and the surgical technique must be based on sound medical judgment by the surgeon that takes into consideration factors such as the circumstances and configuration of the injury.

This document is controlled by TriMed, Inc. When downloaded, printed, and/or copied, this document becomes uncontrolled, and users should always check trimedortho.com for the latest version.